

**ZOLL AED 3®**

**ZOLL®**



THE AED THAT  
**LEADS THE WAY**

# ARE YOU PREPARED TO HELP IN A HEARTBEAT?

It happens in a split second. A person collapses — the victim of sudden cardiac arrest (SCA) — and the clock starts ticking. You call for help and the medical personnel are on their way. But precious seconds are racing by. The victim's best chance for survival is immediate CPR and a rescuer with an automated external defibrillator (AED).



CPR plus early defibrillation can **more than double the rate of survival** from out-of-hospital cardiac arrest.

– 2020 AHA Guidelines<sup>1</sup>

## Designed For Unexpected Heroes™

The latest guidelines from resuscitation councils worldwide are clear: Successful defibrillation must be supported with high-quality CPR. The AED's first heart analysis is critical. If it calls for no shock, only high-quality CPR can lead to a shockable rhythm on the next heart analysis. High-quality CPR provides a struggling heart with oxygenated blood needed to help restore a normal rhythm.

The ZOLL AED 3® defibrillator will guide you through the process of performing high-quality CPR and will, if needed, deliver a potentially lifesaving shock to the heart. With a ZOLL AED 3 and quick action, tragedy could be averted. ZOLL's proprietary Real CPR Help® technology helps guide rescuers in performing high-quality CPR, in accordance with the most recent guidelines for depth and rate of compression. The innovative, intuitive design and enhanced features of the ZOLL AED 3 will give unexpected heroes both the confidence and the knowledge needed to treat SCA.





**“...monitoring of CPR quality is arguably one of the most significant advances in resuscitation practice in the past 20 years and one that should be incorporated into every resuscitation and every professional rescuer program.**


– American Heart Association CPR Quality Consensus Statement<sup>2</sup>

**”**






# EMPOWERING RESCUERS WITH REAL CPR HELP

While only 50% of sudden cardiac arrest victims will initially need a shock, nearly all will require CPR to increase the flow of oxygenated blood to the heart and brain. How can rescuers tell if they are providing high-quality CPR? A ZOLL AED 3 device with Real CPR Help lets rescuers know if they are doing compressions at the most effective rate and depth. The AED's integrated real-time feedback eliminates guesswork; it provides rescuers with visual and audio prompts to "Push Harder" or will let them know when they are doing "Good Compressions," keeping rescuers on track to administer high-quality CPR.




# THE CASE FOR AEDS

Clear, colorful graphic images, together with audio prompts, guide you easily and precisely through the rescue.



## Survival increases with early intervention

Research shows that the probability of survival increases dramatically when CPR is performed and when an AED is applied before an ambulance arrives.<sup>3</sup>



## Where is the AED?

Too often, the answer is, "We don't have one." Sadly, the same research that demonstrated a nearly fivefold increase in survivability (from 5% to 24%) when an AED is used, also showed that an AED is available only 2% of the time.

“

***ZOLL defibrillators equipped with Real CPR Help – providing real-time feedback for depth and rate of chest compressions – combined with training more than doubled the chance of survival from cardiac arrest.***

Bobrow B, et al. Annals of Emergency Medicine<sup>4</sup>

”



# INNOVATIVE DESIGN AND OFFER RESCUERS

## SMART

### Enhanced Real CPR Help

Proven Real CPR Help technology has a full-color display with vivid rescue images, a CPR cycle timer, and a large color bar gauge that shows CPR compression depth. A calm voice prompts you with "Push Harder" or "Good Compressions." The voice will tell you to pause while it completes a heart analysis and will then alert you if a shock is needed or if you should resume CPR until emergency personnel arrive.

### RapidShock™ Analysis

RapidShock analysis enables the industry's shortest rhythm analysis and allows ZOLL AED 3 to deliver a shock, if needed, in as little as 5 seconds. Minimizing pause time allows for more lifesaving CPR, improves CPR quality, and can improve patient outcomes from sudden cardiac arrest.

### Touch screen display

In addition to providing full-color rescue images, text, and real-time CPR feedback, the ZOLL AED 3 screen also includes one-of-a-kind intuitive touch screen capabilities. The touch screen allows quick access to AED information and enables configuration without the need for additional external devices or software.

### Integrated pediatric rescue

The universal design of the CPR Uni-padz™ electrodes gives rescuers a single, versatile solution for both adult and pediatric victims of SCA. To treat a child, use the same set of pads and simply activate child mode. You will be prepared if your victim is seven months old or 70 years old.





# D ENHANCED FEATURES MORE SUPPORT

READY

## WiFi connected

Cloud connectivity can enable automatic reporting of device status, giving you confidence your AED is ready in an emergency. WiFi also provides the ability to quickly access and transmit cardiac arrest event data to medical professionals.

## Long-life consumables

ZOLL's long-lasting batteries and pads increase readiness by reducing the frequency of maintenance.

**CPR Uni-padz:** ZOLL's CPR Uni-padz have an industry-leading shelf life of five years and can be used on both adult and pediatric patients, so there is no need to maintain two separate sets of pads.

**Smart Battery:** The ZOLL AED 3 incorporates a smart battery that can report its power level status on the device's LCD screen or transmit it automatically over WiFi to your AED program management. The battery lasts up to five years once installed and costs less than half the price of competitive smart batteries. This all adds up to savings for your organization and readiness in an emergency.



## HOPE IS IN YOUR HANDS

The ZOLL AED 3 defibrillator represents an invaluable investment in emergency preparedness. ZOLL AEDs give bystanders the confidence to intervene in a case of sudden cardiac arrest and to provide continuous high-quality CPR to increase a victim's chance of survival. With ZOLL AED 3, you'll be ready to help in a heartbeat.

<sup>1</sup>AHA 2020 Guidelines.

<sup>2</sup>AHA CPR Quality Consensus Statement. *Circulation*. 2013;128; pg. 420

<sup>3</sup>Weisfeldt ML, et al. *J Am Coll Cardiol*. 2010;55(16):1713-20.

<sup>4</sup>Bobrow B., et al. *Annals of Emergency Medicine*. July 2013;62(1): 47-56.31

**ZOLL MEDICAL CORPORATION**

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | [zoll.com](http://zoll.com)

Copyright © 2021 ZOLL Medical Corporation. All rights reserved. Designed for Unexpected Heroes, Uni-padz, RapidShock, Real CPR Help, ZOLL AED 3, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.  
MCN PP 2110 0625

For subsidiary addresses and fax numbers, as well as other global locations, please go to [zoll.com/contacts](http://zoll.com/contacts).

**ZOLL**®